

Liberty Nebuliser Pump

Constant use

Product Information



Versatile and Reliable

The Liberty® Nebuliser Pump is the product of 20 years of extensive research and development to provide the patient with the most cost effective, versatile and reliable pump on the market.

The Liberty® Nebuliser Pump has been designed with high performance as a priority. From the air intake filter through to the nebuliser bowl, nothing has been left to chance. The motor and precision pump assembly allows continuous operation without the annoying short-term usage times or cool down periods associated with similarly priced pumps. Use it whenever you want, how often you want and for how long you want.

Having a compact and contemporary design, the Liberty® Nebuliser Pump is a quality product suitable for use with a wide variety of medications recommended by your doctor for the treatment and/or relief of respiratory complaints.

Aerosol therapy is one of the methods of drug delivery for the treatment of chronic and acute pulmonary disease. This form of treatment is increasingly favoured by both family doctors and specialists since it achieves maximum drug efficiency directly in the respiratory tract thus minimising possible sideeffects.

The Liberty® aerosol therapy system incorporates an oil free compressor. Combined with the performance of the RAPIDFLAEM2 jet nebuliser (2.7-3µm) and a full set of accessories including mask, mouth piece and nasal prongs, the Liberty® Nebuliser Pump is a high quality system, making it suitable for basic care and seasonal allergic pathologies.

The Liberty® Nebuliser Pump has a patented anti-vibration system offering a silent and vibration-free operation. The case is moulded from thermoplastic resins, has robust features, a handy carry handle and a large accessory compartment which is hygienically separated from the cord storage compartment.



Features

- Continuously rate for constant use
- Contemporary compact styling
- Super quiet operation
- Compartment for storage of accessories
- Latex free

Useful information

- Continuously rate for constant use
- Contemporary compact styling
- Super quiet operation
- Compartment for storage of accessories

Specifications

Particle Size (Mmad):	2.7 - 3.0 Micron
Power Supply:	230 V ~ 50 HZ / 130 VA
Maximum Pressure:	22 PSI
Flow:	9 LPM
Sound Level:	54 DB
Dimensions:	17 x 31 x 10 cm
Weight:	2.1 KG
Warranty:	2 Years
Approvals:	TGA Artg No. 92850 Class Ila

High Quality Maintenance Free Compressor

High Efficiency Filter



Axis Health

405 – 409 Ferntree Gully Rd,
Mt Waverley, VIC, 3149, Australia

T +61 (3) 9922 0300

F +61 (3) 9544 5392

W www.axishealth.com.au

Code	Description
PNL	Liberty nebuliser pump
PNLF	Spare filter (PK5)
31310	Nebuliser tubing 2 metres
PNLSA	Nebuliser bowl, mask and tubing (set) - Adult
PNLSC	Nebuliser bowl, mask and tubing (set) - Child